

ARG82303
Human HAV antibody ELISA KitPackage: 96 wells
Store at: 4°C

Summary


Product Description	ARG82303 Human HAV antibody ELISA Kit is an Enzyme Immunoassay kit for the qualitative determination of Human HAV antibody in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Human HAV antibody
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	Cut - off
Sample Volume	100 µl
Precision	Intra-Assay CV: less than 20% Inter-Assay CV: less than 25%

Application Instructions

Assay Time	~ 2 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82303 Human HAV antibody ELISA Kit typical data demonstration image

ARG82303 Human HAV antibody ELISA Kit results of a typical data with optical density reading at 450 nm.