

ARG82399 Mouse LYPD3 ELISA Kit

Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG82399-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG82399-002	Standard	2 X 10 ng/vial	4°C
ARG82399-003	Standard/Sample diluent	30 ml (Ready to use)	4°C
ARG82399-004	Antibody conjugate concentrate (100X)	1 vial (100 µl)	4°C
ARG82399-005	Antibody diluent buffer	12 ml (Ready to use)	4°C
ARG82399-006	HRP-Streptavidin concentrate (100X)	1 vial (100 µl)	4°C
ARG82399-007	HRP-Streptavidin diluent buffer	12 ml (Ready to use)	4°C
ARG82399-008	25X Wash buffer	20 ml	4°C
ARG82399-009	TMB substrate	10 ml (Ready to use)	4°C (Protect from light)
ARG82399-010	STOP solution	10 ml (Ready to use)	4°C
ARG82399-011	Plate sealer	4 strips	Room temperature

Summary

Product Description	ARG82399 Mouse LYPD3 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse LYPD3 in serum, plasma (EDTA, heparin, citrate) and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	LYPD3
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	50 pg/ml
Sample Type	Serum, plasma (EDTA, heparin, citrate) and cell culture supernatants.
Standard Range	93.7 - 6000 pg/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: 4.6% Inter-Assay CV: 6.6%

Alternate Names C4.4A; Ly6/PLAUR domain-containing protein 3; GPI-anchored metastasis-associated protein C4.4A homolog; MIG-C4; Matrigel-induced gene C4 protein

Application Instructions

Assay Time ~ 5 hours

Properties

Form 96 well

Storage instruction Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.

Note For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol LYPD3

Gene Full Name LY6/PLAUR domain containing 3


Function Supports cell migration. May be involved in urothelial cell-matrix interactions. May be involved in tumor progression. [UniProt]

Highlight Related products:
[LYPD3 ELISA Kits](#);
New ELISA data calculation tool:
[Simplify the ELISA analysis by GainData](#)

PTM N-glycosylated and O-glycosylated. [UniProt]

Cellular Localization Cell membrane; Lipid-anchor, GPI-anchor. [UniProt]

Images



ARG82399 Mouse LYPD3 ELISA Kit standard curve image

ARG82399 Mouse LYPD3 ELISA Kit results of a typical standard run with optical density reading at 450 nm.