

## ARG82507 Human CNTN5 / Contactin 5 ELISA Kit

Package: 96 wells  
Store at: 4°C

### Component

Cat. No.	Component Name	Package	Temp
ARG82507-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG82507-002	Standard	2 X 10 ng/vial	4°C
ARG82507-003	Standard/Sample diluent	30 ml (Ready to use)	4°C
ARG82507-004	Antibody conjugate concentrate (100X)	1 vial (100 µl)	4°C
ARG82507-005	Antibody diluent buffer	12 ml (Ready to use)	4°C
ARG82507-006	HRP-Streptavidin concentrate (100X)	1 vial (100 µl)	4°C
ARG82507-007	HRP-Streptavidin diluent buffer	12 ml (Ready to use)	4°C
ARG82507-008	25X Wash buffer	20 ml	4°C
ARG82507-009	TMB substrate	10 ml (Ready to use)	4°C (Protect from light)
ARG82507-010	STOP solution	10 ml (Ready to use)	4°C
ARG82507-011	Plate sealer	4 strips	Room temperature

### Summary

Product Description	ARG82507 Human CNTN5 / Contactin 5 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human CNTN5 / Contactin 5 in serum, plasma (EDTA, heparin, citrate) and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	CNTN5 / Contactin 5
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	50 pg/ml
Sample Type	Serum, plasma (EDTA, heparin, citrate) and cell culture supernatants.
Standard Range	62.5 - 4000 pg/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: 4.8% Inter-Assay CV: 6.3%

Alternate Names hNB-2; NB-2; HNB-2s; Contactin-5; Neural recognition molecule NB-2

## Application Instructions

Assay Time ~ 5 hours

## Properties

Form 96 well

Storage instruction Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.

Note For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation

Gene Symbol CNTN5

Gene Full Name contactin 5


Background The protein encoded by this gene is a member of the immunoglobulin superfamily, and contactin family, which mediate cell surface interactions during nervous system development. This protein is a glycosylphosphatidylinositol (GPI)-anchored neuronal membrane protein that functions as a cell adhesion molecule. It may play a role in the formation of axon connections in the developing nervous system. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2011]

Function Contactins mediate cell surface interactions during nervous system development. Has some neurite outgrowth-promoting activity in the cerebral cortical neurons but not in hippocampal neurons. Probably involved in neuronal activity in the auditory system (By similarity). [UniProt]

Highlight Related products:  
[Contactin antibodies](#); [Contactin ELISA Kits](#);  
New ELISA data calculation tool:  
[Simplify the ELISA analysis by GainData](#)

Cellular Localization Cell membrane; Lipid-anchor, GPI-anchor. [UniProt]

## Images



ARG82507 Human CNTN5 / Contactin 5 ELISA Kit standard curve image

ARG82507 Human CNTN5 / Contactin 5 ELISA Kit results of a typical standard run with optical density reading at 450 nm.