

ARG83090
E.coli HCP ELISA KitPackage: 96 wells
Store at: 4°C

Summary

Product Description	E.coli HCP ELISA Kit is designed for the quantitative detection of HCP (host cell protein) content in biopharmaceuticals expressed on E.coli by using a double-antibody sandwich method.
Tested Reactivity	E. coli
Tested Application	ELISA
Target Name	E.coli HCP
Species	E. coli
Conjugation	HRP
Sensitivity	0.5 ng/mL
Standard Range	1 - 243ng/mL

Application Instructions


Assay Time	~2.5 hours
------------	------------

Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Quality control of transgene DNA; Related news: Quality control for engineered T cell therapy
-----------	--



ARG83090 E.coli HCP ELISA Kit standard curve image

ARG83090 E.coli HCP ELISA Kit results of standard run with optical density reading at 450 nm