

ARG83111 HIV1 p24 ELISA Kit (Rapid One-Step)

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83111 HIV1 p24 ELISA Kit (Rapid One-Step) is an Enzyme Immunoassay kit for the quantification of HIV1 p24 in all lentiviral biological samples. It is a rapid One-step 45 minutes protocol with 10ul sample volume.
Tested Reactivity	HIV
Tested Application	ELISA
Target Name	HIV1 p24
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	0.35 ng/mL
Sample Type	All lentiviral biological sample
Standard Range	1.37 - 1000 ng/mL
Sample Volume	10 µl

Application Instructions


Assay Time	~45 minutes
------------	-------------

Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related news: Quality control for engineered T cell therapy
-----------	--



ARG83111 HIV1 p24 ELISA Kit standard curve image

ARG83111 HIV1 p24 ELISA Kit results of standard run with optical density reading at 450 nm