

ARG83386 Adenovirus Hexon ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary


| | |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Description | ARG83386 Adenovirus Hexon ELISA Kit is an Enzyme Immunoassay kit for the quantification of Adenovirus Hexon in Various biological samples types. |
| Tested Reactivity | Other |
| Tested Application | ELISA |
| Target Name | Adenovirus Hexon |
| Conjugation | HRP |
| Conjugation Note | Read at 450 nm. |
| Sensitivity | 0.1 ng/mL |
| Sample Type | Various biological samples types |
| Standard Range | 0.156 - 10 ng/mL |
| Sample Volume | 100 µL |

Application Instructions

| | |
|------------|----------|
| Assay Time | ~5 hours |
|------------|----------|

Properties

| | |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form | 96 well |
| Storage instruction | Store components at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |



ARG83386 Human Adenovirus Hexon ELISA Kit standard curve image

ARG83386 Human Adenovirus Hexon ELISA Kit results of a typical standard run with optical density reading at 450 nm.