

## ARG83563 HbA1c Assay Kit

Package: 100 tests  
Store at: 4°C

### Summary

Product Description	ARG83563 HbA1c Assay Kit can be used to measure HbA1c in Whole blood, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	HbA1c
Conjugation Note	Read at 570 nm
Sensitivity	0.05 µmol/mL
Sample Type	Whole blood, Other biological fluids
Standard Range	0.1 - 5 µmol/mL
Sample Volume	10 µl

### Application Instructions

Assay Time	15 min
------------	--------

### Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related news: <a href="#">Quantifying IR-associated Adipokines with Multiplex ELISA kit;</a>
-----------	---