

**ARG83567**  
**beta HBDH Activity Assay Kit**Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG83567 beta HBDH Activity Assay Kit can be used to measure beta HBDH in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	beta HBDH
Conjugation Note	Read at 450 nm
Sensitivity	1.5 pmol/ml
Sample Type	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	3 - 300 pmol/ml
Sample Volume	10 µl
Alternate Names	EC 1.1.1.30; D-beta-hydroxybutyrate dehydrogenase

### Application Instructions

Assay Time	5 min
------------	-------

### Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.