

ARG83569 Pyruvate Oxidase Activity Assay Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83569 Pyruvate Oxidase Activity Assay Kit can be used to measure Pyruvate Oxidase in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Pyruvate Oxidase
Conjugation Note	Read at 550 nm
Sensitivity	0.025 µmol/ml
Sample Type	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	0.05 - 5 µmol/ml
Sample Volume	10 µl
Alternate Names	EC 1.2.3.3; phosphate-dependent pyruvate oxidase; pyruvic oxidase

Application Instructions

Assay Time	10 min
------------	--------

Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.