

## ARG83580

### arigoQIK® Human IL13 ELISA Development Kit

Package: 1 kit(5 plates), 1 kit (15 plates)  
Store at: 4°C, -20°C

#### Summary

##### Product Description

ARG83580 arigoQIK® Human IL13 ELISA Development Kit, includes Capture antibody, Detection antibody, Standard, and HRP-Streptavidin Solution.  
This ELISA Development Kit is designed for the development of sandwich ELISA to measure Human IL13 in Serum, plasma and cell culture supernatants.

For other reagents required for [arigoQIK® ELISA Development Kit](#), please refer [ARG83524 Integral Reagent Kit \(ELISA Development Kit\)](#)

##### More about arigoQIK®:

- Optimized capture and detection antibody pairs
- Reduced incubation time and wash cycles
- 2-hour quicker than conventional ELISA process
- 5- and 15-plate packages available

##### Tested Reactivity

Hu

##### Tested Application

ELISA

##### Target Name

IL13

##### Conjugation

HRP

##### Conjugation Note

Substrate: TMB and read at 450 nm.

##### Sensitivity

3.9 pg/ml

##### Sample Type

Serum, plasma and cell culture supernatants.

##### Standard Range

7.8 - 500 pg/ml

##### Sample Volume

50 µl

##### Alternate Names

IL13; Interleukin 13; IL-13; P600; Interleukin-13; MGC116786; MGC116788; MGC116789; ALRH; BHR1; Bronchial Hyperresponsiveness-1 (Bronchial Asthma); Allergic Rhinitis; NC30

#### Properties

##### Form

96 well

##### Storage instruction

Store components at 4°C or -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.

##### Note

For laboratory research only, not for drug, diagnostic or other use.

#### Bioinformation

##### Gene Symbol

IL13

##### Gene Full Name

Interleukin 13


##### Background

This gene encodes an immunoregulatory cytokine produced primarily by activated Th2 cells. This cytokine is involved in several stages of B-cell maturation and differentiation. It up-regulates CD23 and MHC class II expression, and promotes IgE isotype switching of B cells. This cytokine down-regulates macrophage activity, thereby inhibits the production of pro-inflammatory cytokines and chemokines.

This cytokine is found to be critical to the pathogenesis of allergen-induced asthma but operates through mechanisms independent of IgE and eosinophils. This gene, IL3, IL5, IL4, and CSF2 form a cytokine gene cluster on chromosome 5q, with this gene particularly close to IL4. [provided by RefSeq, Jul 2008]

Function	Exerts its biological effects through its receptors which comprises the IL4R chain and the IL13RA1 chain, to activate JAK1 and TYK2, leading to the activation of STAT6. Aside from IL13RA1, another receptor IL13RA2 acts as a high affinity decoy for IL13 and mediates internalization and depletion of extracellular IL13. [UniProt]
PTM	Disulfide bond, Glycoprotein. [UniProt]
Cellular Localization	Secreted. [UniProt]

## Images



ARG83580 arigoQIK Human IL13 ELISA Development Kit standard curve image

ARG83580 arigoQIK Human IL13 ELISA Development Kit results of a typical standard run with optical density reading at 450 nm.