

ARG83594
GST ELISA KitPackage: 96 wells
Store at: 4°C, -20°C**Summary**

Product Description	ARG83594 GST ELISA Kit is an Enzyme Immunoassay kit for the quantification of GST (Glutathione S-transferase).
Tested Reactivity	Other
Tested Application	ELISA
Target Name	GST
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	AKR-5185
Sample Type	Cell lysate and Tissue lysate
Standard Range	78.125 - 5000 pg/ml
Sample Volume	100 µL

Properties

Form	Liquid
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.