

ARG83603 SAM + SAH ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C


Summary

Product Description	ARG83603 SAM + SAH ELISA Kit is an Enzyme Immunoassay kit for the quantification of SAM + SAH.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	SAM + SAH
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	SAM: 0.3 nM SAH: 0.1 µM
Sample Type	Plasma, Serum, Cell lysate, Tissue lysate and other biological fluids
Standard Range	SAM: 0.625 - 40 nM SAH: 0.195 - 12.5 µM
Sample Volume	50 µL

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG83603 SAM+ SAH ELISA Kit standard curve image

ARG83603 SAM+ SAH ELISA Kit results of a typical standard run with optical density reading at 450 nm.