

ARG83736 Glucose Assay Kit

Package: 100 tests
Store at: -20°C

Summary

| | |
|---------------------|--|
| Product Description | Glucose Assay Kit is a detection kit for the quantification of Glucose Content in serum, plasma, urine, saliva, milk, culture medium, food and beverages sample. |
| Tested Reactivity | All |
| Tested Application | FuncSt |
| Target Name | Glucose |
| Conjugation Note | Read at 565 nm. |
| Sensitivity | 0.03 mM (0.54 mg/dL) |
| Detection Range | 0.03 - 2 mM (0.54 - 36 mg/dL) |
| Sample Type | Serum, plasma, urine, saliva, milk, culture medium, food and beverages sample. |
| Sample Volume | 20 µl |

Application Instructions

| | |
|------------|--------|
| Assay Time | 30 min |
|------------|--------|

Properties

| | |
|---------------------|---|
| Form | Liquid |
| Storage instruction | Store components at -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |