

ARG83744 Glucagon ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83744 Glucagon ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human, Rat, Mouse, Pig, Bovine, Dog Glucagon in serum and plasma.
Tested Reactivity	Hu, Ms, Rat, Bov, Dog, Pig
Tested Application	ELISA
Specificity	Cross Reactivity: Glucagon 100% Oxyntomodulin 100% Not Cross Reactivity: Glucagon-Like Peptide-1 (GLP-1) 0% Glucagon-Like Peptide-2 (GLP-2) 0% Insulin 0% Somatostatin 0% Amylin 0% CGRP 0% Amylin 0% PPP 0% Pro-Insulin 0%
Target Name	Glucagon
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.16 ng/ml
Sample Type	Serum and Plasma
Standard Range	0.16 - 3.74 ng/ml
Sample Volume	50 µl
Precision	