


ARG83751 Human FGF basic ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83751 Human FGF basic ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human FGF basic in serum, plasma and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	FGF basic
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	7.81 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	15.63 - 1000 pg/ml
Sample Volume	50 µl
Precision	

Images



ARG83751 Human FGF basic ELISA Kit standard curve image

ARG83751 Human FGF basic ELISA Kit results of a typical standard run with optical density reading at 450 nm.