


ARG83753 Human Adiponectin ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83753 Human Adiponectin ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Adiponectin in serum, plasma and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Adiponectin
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	31.25 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	62.5 - 4000 pg/ml
Sample Volume	50 µl
Precision	

Images



ARG83753 Human Adiponectin ELISA Kit standard curve image

ARG83753 Human Adiponectin ELISA Kit results of a typical standard run with optical density reading at 450 nm.