


ARG83754 Human Resistin ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83754 Human Resistin ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Resistin in serum, plasma and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Resistin
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	15.63 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	31.25 - 2000 pg/ml
Sample Volume	50 µl
Precision	

Images



ARG83754 Human Resistin ELISA Kit standard curve image

ARG83754 Human Resistin ELISA Kit results of a typical standard run with optical density reading at 450 nm.