


ARG83761 Rat IL4 ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83761 Rat IL4 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Rat IL4 in serum, plasma and cell culture supernatants.
Tested Reactivity	Rat
Tested Application	ELISA
Target Name	IL4
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	15.63 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	31.25 - 2000 pg/ml
Sample Volume	50 µl
Precision	

Images



ARG83761 Rat IL4 ELISA Kit standard curve image

ARG83761 Rat IL4 ELISA Kit results of a typical standard run with optical density reading at 450 nm.