

ARG83802 General Hyaluronic Acid ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83802 General Hyaluronic Acid ELISA Kit is a detection kit for the quantification of General Hyaluronic Acid in serum, plasma, saliva, urine, cell culture supernatants, cell lysate and tissue lysate samples.
Tested Reactivity	All
Tested Application	ELISA
Conjugation	HRP
Conjugation Note	Read at 450 nm
Sensitivity	7.815 ng/mL
Sample Type	Serum, plasma and other biological fluids
Standard Range	15.63 - 1000 ng/mL
Sample Volume	50 µl

Application Instructions


Assay Time	2.5 hr
------------	--------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related news: TCA intermediate fumarate promotes mitobiogenesis
-----------	--



ARG83802 General Hyaluronic Acid ELISA Kit standard curve image

ARG83802 General Hyaluronic Acid ELISA Kit results of a typical standard run with optical density reading at 450 nm.