

ARG83807 Giardia lamblia ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83807 Giardia lamblia ELISA Kit is a Enzyme Immunoassay kit for the determination of Giardia lamblia in veterinary stool specimens.
Tested Reactivity	Cat, Dog
Tested Application	ELISA
Target Name	Giardia lamblia
Conjugation	HRP
Sample Type	Stool
Standard Range	qualitative
Sample Volume	0.1 g

Application Instructions

Assay Time	1 h (37°C), 30, 15 min (RT)
------------	-----------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Microbiology and Infectious Disease kit
---------------	---