

ARG90006 WB Incubation Box

Package: 1 EA

Summary

Product Description	ARG90006 WB Incubation Box is a useful tool for antibody incubating during WB process. The incubation box including four chambers and each chamber requires as little as 2 ml of solution for WB incubation. (SIZE: 12 X 9 X 2 cm)
Tested Application	WB

Images



ARG90006 Incubation Box image

ARG90006 Incubation Box with 4 chambers and a sealable lid.



ARG90006 Incubation Box image

ARG90006 Incubation Box is a convince incubation box for antibody incubation during western blot process. Using this box, scientist could handle 4 antibodies/mini-membranes at once and it requests as little as 2 ml for antibody incubation.